

Needlestick prevention that meets your needs



virtually eliminates exposure ◀

single-handed, passive activation ◀

hands remain behind the sharp ◀

integrated safety mechanism ◀

prevents recapping and needle removal ◀

remains safe through disposal ◀

easy to use, requires minimal training ◀

safe and effective for patients ◀

All syringes are packaged 100 syringes per box

Catalog Number Description

1cc VanishPoint® Tuberculin syringe

10131 1cc VanishPoint syringe 27G x 1/2"
10151 1cc VanishPoint syringe 25G x 5/8"
10161 1cc VanishPoint syringe 25G x 1"

1cc VanishPoint® U-100 Insulin syringe

10211 1cc VanishPoint syringe 29G x 1/2"

3cc VanishPoint® syringe

10301 3cc VanishPoint syringe 25G x 5/8"
10391 3cc VanishPoint syringe 25G x 1"
10311 3cc VanishPoint syringe 23G x 1"
10331 3cc VanishPoint syringe 22G x 1"
10341 3cc VanishPoint syringe 22G x 1 1/2"
10351 3cc VanishPoint syringe 21G x 1"
10361 3cc VanishPoint syringe 21G x 1 1/2"
10371 3cc VanishPoint syringe 20G x 1"
10381 3cc VanishPoint syringe 20G x 1 1/2"

5cc VanishPoint® syringe

10531 5cc VanishPoint syringe 22G x 1"
10541 5cc VanishPoint syringe 22G x 1 1/2"
10551 5cc VanishPoint syringe 21G x 1"
10561 5cc VanishPoint syringe 21G x 1 1/2"
10571 5cc VanishPoint syringe 20G x 1"
10581 5cc VanishPoint syringe 20G x 1 1/2"

10cc VanishPoint® syringe

11031 10cc VanishPoint syringe 22G x 1"
11041 10cc VanishPoint syringe 22G x 1 1/2"
11051 10cc VanishPoint syringe 21G x 1"
11061 10cc VanishPoint syringe 21G x 1 1/2"
11071 10cc VanishPoint syringe 20G x 1"
11081 10cc VanishPoint syringe 20G x 1 1/2"

VanishPoint® blood collection products

22701 VanishPoint blood collection tube holder
250 per box
22711 VanishPoint small diameter tube adapter
25 per box



**RETRACTABLE
TECHNOLOGIES, INC.**

www.vanishpoint.com

CE
0044

©2003 Retractable Technologies,
Inc. All rights reserved.M005 REV. 2

622 South Mill Street
Lewisville, Texas
75057-4632

Tel: (972) 221-6644
Fax: (972) 221-9786
Toll Free: 1-888-703-1010

UNCOMPROMISING
SAFETY FOR CLINICIANS
AND PATIENTS



VANISHpoint™
SYRINGES

The New Standard for Safety™

VANISHpoint™

1

Use standard procedure
for drawing up medication.



2

Prepare and give injection
using aseptic technique
according to institutional policy.



3

While the needle is still in
the patient, fully depress the
plunger to activate retraction.



4

Dispose of VanishPoint®
syringe in an appropriate
sharps container.

